



**KERALA MEDICAL SERVICES CORPORATION LTD**

(A Government Of Kerala Undertaking)

Thycaud.P.O., Thiruvananthapuram-14.

Tele Fax No. 0471-2945647

0471-2945646, 0471-2945600

E-mail id : [etenders\\_ep@kmscl.kerala.gov.in](mailto:etenders_ep@kmscl.kerala.gov.in)

CIN : U24233KL2007SGC021616

GSTIN: 32AADCK4029M1ZK

**Equipment Division**

[www.kmscl.kerala.gov.in](http://www.kmscl.kerala.gov.in)

**No: KMSCL/EP/4058/2022**

**Dated: 29.07.2022**

**NOTICE**

**Sub:-** Tender No.T437/2021 - Cancellation Notice –Low Temperature Plasma Sterlizer 100 Litre- reg.

**Ref:-** 1. Tender Inviting Notice published dated 25.01.2022

2. Technical bid opening held on 28.02.2021 & 11.03.2022 @ 11.00 A M

3. Eligible list published on 25.05.2022

4. Reports of demonstration/ document evaluation of Low Temperature Plasma Sterlizer 100 Litre received on 21.07.2022

I. The reports of demonstration/ technical document evaluation of : Low Temperature Plasma Sterlizer 100 Litre as per reference no. 4th cited above was submitted by the technical committee appointed for the purpose. The details of the tenderers qualified in the technical evaluation are tabulated as follows:

Sl. No	Equipment	Name of the Selected Bidder	Selected Models
1.	:Low Temperature Plasma Sterlizer 100 Litre	M/s. Vibrant healthcare solutions	STERRAD 100S/Advanced sterilization

II. The other tenderers participated in the tender are rejected due to the reasons mentioned against their name below;

1. M/s. CL Micromed Pvt. Ltd.

Rejected Model: Lowtem Crystal 120M/LowTem Co

Reason for rejection: Sample not Produced

2. M/s. Genist technocracy vt ltd

Rejected Model: GTPS-001/Genist Technocracy

Reason for rejection: Sample not produced

3. M/s. Kerala Surgical Equipment Co

Rejected Model: Eiligplaz/Sterile Safequip

Reason for rejection;

Clause No	Technical Specification	Deviations Noted
4	Machine should operate with standard 3 phase supply.	Single Phase Only
5	The sterlant should be cassette based H2O2 injection system with approximate 2ml per injection. Hydrogen peroxide >55% sterlant should be in cassette with leak proof indicator and H2O2 should be transported by AIR mode for quick supplies. The rate of 1 cassette should be offered in BOQ and for L1 evaluation rate of 100 nos shall be taken by the BOQ.	Bottle Type
6	Should have the capacity to carry more than 40kg of load per cycle.	Below 25 Kg Only

III Due to Administrative reason the tender for **KMSCL/EP/T437/1515B/2022: Low Temperature Plasma Sterilizer 100 Litre** is cancelled.

**Dr. S S Joy**  
General Manager